B-Med

B-HIP IMPLANT SYSTEM

B-Med S.r.l.

Via Bazzanese 32/7
CAP 40033
Casalecchio di Reno (BO)
P.IVA 03800221206
Mail: info@b-med.eu
PEC: info@pec.b-med.eu
Web: www.b-med.eu

I

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GB

B-HIP IMPLANT SYSTEM – RECOMMENDATIONS AND PRECAUTIONS

General remarks

It is advisable that the patient is fully informed by the surgeon of the reasons recommending a hip joint prosthesis implant, of its advantages as well as of its contra-indications and possible post-op complications.

The patient should also be informed of the type of surgery that will be practised and of the expected functional and aesthetic outcome.

Before proceeding with the operation, a final and thorough evaluation must be made taking the pros and cons of this type of surgery into consideration.

In particular, attention must be given to the patient's overall clinical condition and the operation must be thoroughly planned beforehand and performed in the best possible manner to ensure a successful outcome.

The following points should be kept in mind:

- an artificial joint is not as functionally efficient as a natural one, but in any case, the post-surgical outcome can be expected to represent an improvement compared to the pre-surgical condition;
- wear, infection, prolonged stress and strain are all factors that may lead an artificial joint to become mobilised and a follow-up operation may be required;
- in addition to other factors, the patient's biological age should be closely considered in deciding in favour of a cemented or non-cemented implant;
- the prosthesis must not be submitted to excessive stress and strain by overweight and/or physically over-active patients;
- the patient must be made to understand and appreciate that the best outcome also depends on taking certain personal precautions and complying with medical prescriptions;
- if intolerance to the presence of foreign bodies is suspected, appropriate tests should be run to make sure that the materials to be used in the operation are well tolerated.

The materials used for manufacturing joint prosthesis components are clearly specified on the package label together with the reference standard, where applicable. The labels inside the package are intended to ensure product traceability and should therefore be affixed to the patient's medical record and/or to the operation report.

The operation must be carefully planned beforehand in order to determine the size of the implant components and their positioning in the bone. A checklist of all the components required for the operation and their availability in the operating theatre should be made and the working efficiency of implant instruments should be ascertained. It should also be borne in mind that apart from the intrinsic quality of the material a prosthesis' working life also depends on a number of biological and biomechanical factors.

ADVISE THE PATIENTS THAT, BEFORE UNDERGOING ANY IMAGING PROCEDURES, THEY NEED TO INFORM THEIR RADIOLOGIST THAT THEY HAVE BEEN IMPLANTED WITH AN ORTHOPAEDIC PROSTHESIS

Indications

- Extensive primary and secondary destruction of the joint to the extent that the functional efficiency of the locomotive apparatus is reduced;
- Severe pathological condition affecting the articulation caused by degenerative and rheumatoid arthritis;
- Joint fracture or bone necrosis;
- Post-surgical conditions after previous operations with or without consequent use of a prosthesis.

Contraindications

- On-going inflammatory process in the periarticular region;
- Severe loss of bone tissue such as to inhibit a primary stabilisation of the prosthesis;
- Degenerative changes in the patient's neurological condition;
- Severe instability in the ligament area that cannot be remedied;
- Foreseeable causes of fatigue of the implanted joint due to obesity or excessive physical activity;
- Severe osteoporosis;
- Bone cancer in the implant anchoring area;
- Alcohol and drug abuse;
- Allergy to the materials employed;
- Lack of collaboration by the patient.

Relative Contraindications

- Adiposity
- Lacking or foreseeable not assured compliance
- Foreseeable overload/overstressing of the joint prosthesis
- Osteoporosis

Complications:

- hematomas in the region of the operation;
- late onset of acute infections in the region of the operation;
- momentary or persistent functional alterations in the nerves of the anatomical area concerned;
- venous thrombosis, pulmonary embolism, heart failure;
- change in position and/or loosening of the prosthesis;
- joint dislocation;
- shortening or lengthening of the limb concerned;
- pathological bone fracture caused by changes in load;
- allergic reactions or metallosis in the peripheral region of the implant;
- periarticular ossification.

Product handling

Sterile products must be stored in their protective packaging at room temperature.

A careful visual inspection of the device should be made before use to check against any damage. Ceramic products believed to have suffered an impact (accidental fall, for instance) must never be used.

Sterilisation

Utmost care is taken by the manufacturer to ensure the maximum safety of use and the best surgical outcome of its implants.

It is likewise essential that operators be just as diligent in handling the prostheses and in implanting them so as to guarantee a successful outcome, taking care to also comply with the recommendations for ensuring the necessary aseptic conditions.

It is fundamental to respect the following rules:

Sterile products are delivered in packages marked STERILE on the label.

The products may have been sterilised by gamma/beta radiation or ethylene oxide (EO); the sterilisation method used is specified on the label.

Before using sterile material always remember to check:

- the sterilisation expiry date (year and month) marked on the label;

(Never use products with an expired sterilisation date.)

- that the package is intact;

(Never use a product from a ruptured package.)

Special instructions for use

- Implant surfaces must not be damaged. It is strictly forbidden to use damaged implants, or implants that may have accidentally fallen or that have already been used or removed during the operation. Coated implants that are stained with blood or smeared by secretory matter must never be re-used.
- Never exert excess force and avoid using metal hammers to fit bioceramic heads into place. Simply tap gently with a plastic hammer. Make sure that outer and inner cone surfaces are perfectly clean and dry.
- If after fixing to the stem cone, the spherical bioceramic head needs to be removed, it can no longer be used.

Fixing of the head risks being inadequate if a damaged cone is used as abrasion may be caused between the prosthetic stem cone and the sphere.

What to do in case of a broken spherical ceramic head:

Proceed as follows, if a spherical ceramic head breaks: carefully remove all ceramic fragments. The polyethylene cotyloid cavity will have to be replaced, even if in good condition. In fact, minute ceramic fragments undetectable to the naked eye may increase wear of the various parts.

Other general instructions:

Careful pre-surgical planning and the adoption of as gentle a surgical technique as possible are a must to prevent the risk of bone fractures that may occur particularly when scraping the femoral bone with rasps when preparing it for the implant or when inserting the implant.

If cement is used always take care to remove any cement that may with time detach itself.

Coupling

Stem-head coupling

Coupling must be made solely between stems and heads both produced by B-MED and having the same type of taper configuration.

— B-MED stems with 12/14 cone ↔ B-MED heads with 12/14 cone

Coupling of femoral cup, insert with head

B-MED Polyethylene femoral inserts can be coupled with B-MED bioceramic and metal joint heads of matching diameter (example: a cup/insert with a spherical cavity of 28 mm with a 28-mm external diameter head) made according to ISO 7206/2 Standards.

Dual Mobility Inserts couple exclusively with Dual Mobility B-MED liners.

The B-MED MS Cups mate with the bone screws of B-MED System.


In the presence of abrasive powders, it is recommended to flush with saline the surgical bed and aspirate solution with suction unit. The utmost care is required in the head-insert coupling so as to avoid the arise of any possible pathology due to the wear and tear of materials. It also recommended to inform the patient about the possible risks of intoxication relating to the revision surgery. Never revise an articular solution after rupture of a ceramic component with a metal-polyethylene articular solution: in such a case adopt a ceramic-polyethylene coupling solution. Never replace a ceramic head with another ceramic head but with special ceramic heads properly designed for revision.

Modular acetabular parts

Coupling instructions are given on a leaflet inside the package and on the package label.

B-MED does not assume any responsibility for any consequences deriving from any couplings other than those expressly indicated by the manufacturer.

Symbols

If in doubt as to any of the instructions and recommendations set out in the present document, please contact the manufacturer who is always available for any clarifications.

Risks associated with re-use of single use devices:

B-MED does not assume any responsibility for any consequences deriving from any re-use of SINGLE USE devices. The manufacturer alerts the end-users of the following risks associated with re-use of SINGLE USE devices:

- Cross-infection transmission;
- Mechanical / Dimensional features damages as a consequence of first use, no warranty of mechanical / dimensional conformity maintenance;

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