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Product handling

Sterile products must be stored in their protective packaging at room temperature. A careful visual inspection of the device should be made before use to check against any damage. Ceramic products believed to have suffered an impact (accidental fall, for instance) must never be used.

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Sterilization

It is taken by the manufacturer to ensure the maximum safety of use and the best surgical outcome of its implants. The sterilization must be just as diligent in handling the prostheses and in implanting them so as to guarantee a successful outcome, taking care to also comply with the recommendations for ensuring the necessary aseptic conditions. It is fundamental to respect the following rules:

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Special instructions for use

- Implant surfaces must not be damaged. It is strictly forbidden to use damaged implants, or implants that may have accidentally fallen or that have already been used or removed during the operation. Coated implants that are stained with blood or smeared by secretory matter must never be re-used.
- Never exert excess force and avoid using metal hammers to fit bioceramic heads into place. Simply tap gently with a plastic hammer. Make sure that outer and inner cone surfaces are perfectly clean and dry.
- If after fixing to the stem cone, the spherical bioceramic head needs to be removed, it can no longer be used.

Fixing of the head risks being inadequate if a damaged cone is used as abrasion may be caused between the prosthesis and the sphere.

What to do in case of a broken spherical ceramic head:

- Proceed as follows, if a ceramic head breaks: carefully remove all ceramic fragments. The polyethylene cotyloid cavity will have to be replaced, even if in good condition. In fact, minute ceramic fragments undetectable to the naked eye may increase wear of the various parts.

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EN Pre-Dicate-HIP IMPLANT SYSTEM – RECOMMENDATIONS AND PRECAUTIONS

General remarks

It is advisable that the patient is fully informed by the surgeon of the reasons supporting a hip joint prosthesis implant, its advantages as well as of its contraindications and possible post-op complications.

Coupling of femoral cup/cotyloid insert with head

LINCOTEK BOLOGNA bicromatica femoral inserts can be coupled with LINCOTEK BOLOGNA bioceramic and metal joint heads of matching diameter (example: a cup/insert with a spherical cavity of 28 mm with a 28-mm external diameter head) made according to ISO 7206/2 Standards.

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Symbols

LOT	EXPIRY DATE	DO NOT REUSE
▲ See instructions for use	REF Reference number	STERILE EO
Product sterilised by ethylene oxide		Product sterilised by radiation

Indicazioni

- Extensive primary and secondary destruction of the joint to the extent that the functional efficiency of the locomotive apparatus is reduced;
- Severe pathological condition affecting the articulation caused by degenerative and metabolic arthritis;
- Severe instability in the ligament area that cannot be remedied;
- Fracture of the bone or cartilage;
- Post-surgical conditions after previous operations with or without consequent use of a prosthesis;

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Indicazioni

- Extensive primary and secondary destruction of the joint to the extent that the functional efficiency of the locomotive apparatus is reduced;
- Severe pathological condition affecting the articulation caused by degenerative and metabolic arthritis;
- Severe instability in the ligament area that cannot be remedied;

LINCOTEK BOLOGNA, does not assume any responsibility for any consequences deriving from any reuse of SINGLE USE devices. The manufacturer alerts the end-user of the following risks associated with re-use of SINGLE USE devices:

- Cross-Infection transmission;
- Mechanical / Dimensional features damages as a consequence of first use, no warranty of mechanical / dimensional conformity maintenance;

Contraindicationen

- On-going inflammatory process in the periaricular region;
- Severe loss of bone tissue such as to inhibit a primary stabilisation of the prosthesis;
- Degenerative changes in the patient's neurological condition;
- Severe instability in the ligament area that cannot be remedied;

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Relative Contraindications

- Adiposity;
- Lacking or foreseeable not assured compliance;
- Foreseeable overload/overstressing of the joint prosthesis;
- Osteoporosis;

Relative Contraindicationen

- Adiposity;
- Lacking or foreseeable not assured compliance;
- Foreseeable overload/overstressing of the joint prosthesis;
- Osteoporosis;

Arthroskopische Operationsmethode:

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