G1

G3

DE

Röntgenpositiver chirurgischer Zement

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Main characteristics

Intended use

G1 and **G3** are radiopaque bone cements for surgical use to perform arthroplasty procedures, such as hip replacement, knee replacement, ankle replacement, shoulder replacement and other joint replacements, specifically formulated to allow fixation of prosthetic devices on the living bone.

G1 and **G3** acrylic-based bone cements demonstrate a short setting time and fast hardening, extremely high mechanical performance, high intrusion, low heat generation rate (versus ISO 5833:2002 – “Implant for surgery-acrylic resin cements”).

G1 standard viscosity allows manual spatula application and comes in a 20 g and 40 g powder package (**G1 20** and **G1 40**).

G3 low viscosity allows syringe application and comes in a 20, 34.4, 60 g powder package (**G3 20**, **G3 40** and **G3 60**).

G1 and **G3** bone cements must be prepared immediately before application, using the two sterile components.

Preparation, handling and application of **G1** or **G3** bone cements must be performed only by qualified healthcare professionals, specifically trained in the procedure and under the direct supervision of the physician responsible for the procedure.

The manufacturer declines any liability in case of the use of **G1** or **G3** bone cements not strictly corresponding to the specified intended use, or the use of **G1** or **G3** bone cements by personnel not adequately qualified and trained.

Composition

Powder components
PMMA (Polymethyl Methacrylate)
BPO (Benzoyl Peroxide)
Barium Sulfate
Liquid components
MMA (Methyl Methacrylate)
DMPT (N,N dimethyl-p-toluidine)
MEHQ

Content

	G1		G3		
	20	40	20	40	60
Powder weight (g)	20	40	20	34.4	60
Liquid volume (ml)	9.1	18.3	11.6	20.0	35.0 ¹

Packaging and sterility

The manufacturing and packaging process of **G1** and **G3** is performed under strict quality procedures in a controlled environment, which conforms to the most stringent international standards.

The liquid component is sterilized by filtration, the powder component is sterilized by ethylene oxide. The sterile liquid

is contained in an amber glass vial packed in a blister pack in sterile conditions and subsequently sterilized by ethylene oxide. The powder is packed in two bags in sterile conditions. The internal bag in medical-grade paper and polyethylene containing the powder component is inserted in another bag in Tyvek and polyethylene and both are sterilized by ethylene oxide. The two bags are packed in a non-sterile protective aluminium wrapping.

Before using **G1** or **G3**, carefully check the integrity of the packaging. If the package is incomplete, damaged or unsealed, the cement cannot be used and must be discarded. After opening the package, it is mandatory, and the responsibility of the operator, to use an aseptic handling technique. Any error in handling and during transfer into the sterile field might affect bone cement sterility, the sterility of the surgical intervention and imply a risk of severe complications for the patient, such as infections and sepsis.

Clinical indications

Bone cements **G1** and **G3** are intended to be used during surgical interventions caused by:

- Arthroplasty revision
- Joint replacement after a traumatic event
- Osteoarthritis
- Osteoporosis hangover

Contraindications

Do not use **G1** and **G3** when a patient is known to be hypersensitive to the constituents of bone cement or to the contrast medium (barium sulfate).

Absolute contraindications are: local or systemic infections not completely resolved; application where there is no possibility of bone regeneration.

Side effects

- It is commonly known that all bone cements may potentially and rarely directly or indirectly cause adverse patient reactions affecting the cardiovascular system. Hypotensive reactions have occurred between 10 and 165 seconds following the application of bone cement; they have lasted from 30 seconds to 5 or more minutes. Some have progressed to cardiac arrest. Patients should be monitored carefully for any change in blood pressure during and immediately following the application of bone cement.
- Rare severe adverse events, associated with the use of bone cements include myocardial infarction, cardiac arrest, cerebrovascular accident, pulmonary embolism, thrombophlebitis, haemorrhage, haematoma, infection of surgical wounds.
- Use of bone cement during pregnancy and breastfeeding: no contraindication known, but use of bone cement should be avoided.
- The patient must be made aware and accept potential risks and complications connected with the use of bone cement. The patient needs to be

¹ In two 17.5 ml vials to be used at the same time.

informed about the measures that could be taken to reduce possible negative effects.

Familiarization with G1 and G3 bone cements

Before using **G1** or **G3**, the surgeon should become familiar with their properties, handling, application, and with the handling of other devices which come into contact with the bone cement, or are going to be used. The surgeon shall perform simulations of the procedure following correct medical practice and the instructions of the manufacturers of involved devices. It is also recommended for surgeons to practice mixing, handling and application prior to use. The characteristics of the bone cement and the timing of the different phases depend on the ambient temperature, humidity and mixing technique. The surgeon shall simulate the entire process of mixing, handling and application in similar environmental conditions to the operating room.

Warnings when using G1 or G3 bone cements

- Monitor the patient carefully for any change in blood pressure during and immediately following the application of the bone cement. Adverse patient reactions affecting the cardiovascular system have been associated with the use of bone cements. Hypotensive reactions have occurred between 10 and 165 seconds following the application of bone cement; they have lasted from 30 seconds to 5 or more minutes. Some have progressed to cardiac arrest. Patients should be monitored carefully for any change in blood pressure during and immediately following the application of bone cement.
- Rare severe adverse events, associated with the use of bone cements include myocardial infarction, cardiac arrest, cerebrovascular accident and pulmonary embolism.
- In the case of cardiovascular or pulmonary complications it is necessary to control and re-establish an appropriate level of volemia.
- Be careful during the bone cement mixing phase and strictly follow the instructions.
- Avoid direct skin or eye contact with the liquid component of dough and reduce as far as possible exposure to monomer vapors, which may cause irritation of the airways, eyes and, in rare cases, affect the liver.
- Ventilate the room to eliminate monomer vapors. The liquid component is volatile and flammable. In the presence of monomer vapors do not use electrocautery instruments or other high temperature sources.
- Do not use latex gloves or other latex devices. The liquid component is a lipid solvent which can cause glove perforation and may damage exposed tissues. PVP gloves (three layers: polyethylene, vinyl copolymer, polyethylene) or Viton®-butyl gloves give adequate protection for a long period. If surgical synthetic rubber gloves are

going to be used, it is advisable to wear a second pair of gloves on top, suitable for bone cement handling.

- Operators wearing contact lenses shall not mix bone cement or be exposed to monomer vapors.
- The polymerization of bone cement is an exothermic reaction which finishes only when the bone cement becomes hard inside the body. Respect the waiting time before application (fig. 1) and consider that generated heat may damage bone tissue or other tissues in contact with bone cement.
- **Do not add any other material to G1 or G3 bone cements.** The addition of non-approved ingredients (powder, water solutions,...) severely compromises physical and chemical characteristics of bone cement both in the preparation phase and after the implant.
- Use always the entire content of one powder bag mixed with the entire content of one single vial (two vials in the case of **G3 60**). It is not allowed to mix more than one powder dose at a time.
- The operator should be instructed on bone cement application techniques, and should respect them strictly. The implant bone cavity should be carefully cleaned, irrigated and dried before bone cement application to avoid contamination of the bone cement through blood or bone marrow. The graft of a foreign body in contact with biological tissues increases the risk of infection in the post-surgical phase.
- Single use device. This device cannot be re-used.
- The safety and effectiveness of G1 and G3 bone cement in pregnant women or children has not yet been established. Bone cement must not be used during the first trimester of pregnancy. During all other stages of pregnancy, the product must only be used in the case of life-threatening circumstances.
- When treating children, G1 and G3 bone cement can be used to save limbs only in the event that other treatment procedures will not have a good probability success.

Application instructions

When using **G1** or **G3**, respect normal application rules for bone cements to limit unpredicted effects. The bone cement implant needs to be performed by a specialist physician, specifically trained.

Preliminary work process

To the prepare bone cement, the following is necessary:

- a sterile working surface;
- a sterile bowl made of ceramic, stainless steel, polypropylene or other material specifically approved for contact with bone cement dough;
- a sterile spoon or spatula made of ceramic, stainless steel, polypropylene or other material specifically approved for contact with bone cement dough.

Only after checking the integrity of the entire packaging can the external non-sterile aluminium laminated envelope be opened by an assistant who will then open only the contained outer pouch, which is the sterility barrier, and then the internal pouch, containing the powder, can be laid on the sterile working surface using an aseptic technique. Similarly the assistant will open the blister containing the vial and then place the glass vial in the sterile field using an aseptic technique.

Any error in handling and during transfer into the sterile field might affect bone cement sterility, the sterility of the surgical intervention and imply a risk of severe complications for the patient, such as infections and sepsis.

A dose is prepared by pouring the entire content of the liquid component vial (two vials in case of **G3 60**) into a bowl containing the entire quantity of powder included in one bag.

Preparation and application

The preparation and application of **G1** or **G3** bone cements take place in four subsequent phases:

- I. mixing;
- II. waiting (in this phase the bone cement, when touched, sticks to the gloves);
- III. application (in this phase the bone cement has the aspect of a fluid dough, it does not stick and can be handled, injected by a syringe and a cannula);
- IV. setting (in this phase the bone cement should have already been implanted and hardening completed. Non-implanted bone cement has in this phase becomes too viscous, it cannot be used anymore and must be discarded).

The duration of phases II to IV depends on the ambient temperature and humidity. A higher temperature accelerates hardening, while a lower temperature slows it down (see fig. 1).

Prepare the mix choosing the appropriate type of bone cement and the dose which is necessary for the kind of surgical intervention and the volume of the cavity to be filled. Always open the phial at a safe distance from the container (risk of shards of glass). All the liquid must be poured into the container first.

Carefully open the bag and pour the entire contents into the container.

The powder must always be added to the liquid.

Mix the entire quantity of powder and liquid in the container well, until you obtain an even mixture that is free of air bubbles.

Mix by hand inside the sterile bowl using a sterile spoon or spatula conforming to materials' specifications.

If you use a mixing or application device for surgical cements, follow the manufacturer's instructions.

Use a regular stirring action, not too fast, and continue for one minute. Do not exceed the mixing time.

Viscosity increases progressively as a consequence of the polymerization reaction, during phases II to IV.

Keep waiting until the completion of phase II (waiting phase) and then proceed with the application (phase III, the cement has become dough).

Dose

Choose the dose considering the intervention to be performed. As a preventive measure, it is advisable to keep at hand at least a second package of **G1** or **G3** bone cement to be able to prepare an extra quantity of bone cement, if required during the procedure.

The maximum recommended dose of cement is as follows:

Pack	No. of packs
G1 20	6
G1 40	3
G3 20	6
G3 40	3
G3 60	2

Other uses

G1 and **G3** are not permitted for any use, other than the specified intended use.

Storage

G1 and **G3** bone cements must be stored in the sealed original package, inside a dry, clean storage room, at a storage temperature not below 5°C and not above 25°C.

Expiry date, sterility, disposal

The expiry date is written on the external box and on the secondary outer blister labels.

G1 and **G3** cannot be used after the expiry date. **G1** and **G3** cannot be used if the package is unsealed or the secondary blister is not present or damaged.

Opened or damaged **G1** and **G3** packages must be disposed of with all their content. **G1** and **G3** are sterilized by ethylene oxide, aseptic filtration and radiation and cannot be re-sterilized. A yellowish colour of the powder or of the liquid is not normal, in this case the cement cannot be used and must be discarded. Dispose of the remaining **G1** or **G3** bone cements and the content of partially utilized, expired, un-usable or damaged packages, following rules and procedures applicable to this kind of hospital waste.

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G1 40 Radiopaque bone cements

Device Description

G1 40 bone cement is a self-curing, radiopaque, polymethylacrylate (PMMA) based cement which is used for securing a metal or polymeric prosthesis to living bone in arthroplasty procedures. The hardened bone cement secures the fixation of the grafted artificial joint improving the transfer of forces at the interface of the implant-bone.

The bone cement is supplied as a two-component system, consisting of separate, sterile liquid and powder components, which are mixed together at the point of use to produce the cement.

The liquid component is sterilized by membrane filtration and aseptically filled into a sterile glass ampoule. The ampoule is contained within a sealed blister pack, which is sterilized using ethylene oxide. The powder component is contained in a Tyvek/PE pouch, within a peelable pouch and is sterilized by ethylene oxide. The sterile powder component is supplied within an outer, protective, non-sterile, laminated foil pouch.

G1 40 bone cement is a standard viscosity cement, primarily intended for manual application.

Pack size	
Powder weight (g)	40,00
Liquid volume (ml)	18,30

Composition

Composition of the powder component (% w/w)	
Polymethyl Methacrylate	88,8
Benzoyl Peroxide	1,2
Barium Sulfate	10,0

Composition of the liquid component (% w/w)	
Methyl Metacrylate	97,50
N,N dimethyl-p-toluidine	2,5
Methy Ether of Hydroquinone	50 ppm

Intended use

G1 40 radiopaque bone cements is intended for use in arthroplasty procedures of the hip, knee, ankle, shoulder and other joints for the fixation of polymer or metallic prosthetic implants to living bone.

Contraindications

The use of **G1 40** is contraindicated:

- when a patient is known to be hypersensitive to the constituents of bone cement or to the contrast medium (barium sulfate);
- when the local or systemic infections not completely resolved;
- where there is no possibility of bone regeneration;
- during pregnancy or breast-feeding: The safety and effectiveness of the radiopaque bone cements in pregnant women or in children has not yet been established.

Precautions/Warnings

Bone cement preparation:

- **G1 40** bone cement is supplied sterile for single use only. Do not re-use. Re-sterilization of any components of the cement must not be attempted.
- Ensure the inner packages and components are undamaged. Powder should be consistent (no agglomerations) and not yellow or brown in color. The contents within the vial should appear as a low viscosity liquid. Do not use the liquid monomer if it shows any sign of thickening or premature polymerization. Always check the condition of the liquid monomer before performing the procedure. If the powder has a yellowish or brownish color or if the liquid is syrupy, do not use product. This indicates the product has not been stored properly.
- A dose is prepared by mixing the entire contents of the bag of cement powder with all the monomer liquid of an ampoule. The quantity of cement dough required depends on the specific surgical intervention and on the technique being used. At least one additional dose of **G1 40** should be available before commencing the operation.
- The protective outer foil pouch and the outer peelable pouch are non-sterile and must not be transferred into the sterile field
- The surgeon should, by specific training and experience, be thoroughly familiar with the properties, handling characteristics and application of bone cements. The user is advised to practice the entire procedure of mixing, handling and introducing **G1 40** before using it for the first time. Detailed knowledge is necessary even if mixing systems and syringes are being used for application of the cement. Because the handling and curing characteristics of these cements vary with temperature and mixing technique, they are best determined by the surgeon's actual experience.
- Handling characteristics and setting time of bone cements are affected by:
 - o temperature
 - o humidity
 - o vacuum mixing
- Refer to clinical timing and usage charts at the end of this leaflet for effect of temperature on handling and setting times.

- Manual handling of the cement and body temperature will reduce the final setting time.
- Variations in humidity will affect the cement handling characteristics and setting time.
- Variations in setting time over the cement's shelf life can be minimized by storing the cement under the recommended conditions
- The surgeon should be familiar with the effects of a particular mixing system on handling characteristics and setting time of the cement before use. Vacuum mixing can noticeably accelerate the setting time of the product.

Handling of the liquid monomer:

- As the liquid monomer is highly volatile and flammable, the operating room should be adequately ventilated to eliminate as much monomer vapor as possible. Ignition of monomer fumes caused by use of electrocautery devices in surgical sites near freshly implanted bone cements has been reported.
- Caution should be exercised during the mixing of the two components to prevent excessive exposure to the concentrated vapors of the monomer, which may produce: irritation of the respiratory tract, eyes, and possibly the liver; contact dermatitis. If the liquid component comes into contact with the eyes, wash with copious amounts of water.
- Concentrated vapors of the liquid component may have an adverse reaction with contact lenses. Since soft contact lenses are permeable to liquids and gases, they should be worn in the operating theater if methyl methacrylate is being used.
- Methyl methacrylate has been demonstrated to cause hypersensitivity in susceptible persons, which may result in an anaphylactic response.
- The liquid component of bone cement is a powerful lipid solvent. This liquid component should not be allowed to come into contact with rubber or latex surgical gloves. Wearing of a second pair of gloves and strict adherence to the mixing instructions may diminish the possibility of hypersensitivity reactions. Gloves made of PVP (polyethylene, ethylene vinylalcohol copolymer, polyethylene) and Viton ®/butyl gloves have proved to provide good protection over a lengthy period. For safety's sake it is recommended that two pairs of gloves be worn over one another, e.g. one polyethylene surgical glove over an inner pair of standard latex surgical gloves. The use of latex or polystyrene-butadiene gloves on their own is inadequate. Please make enquiries with your supplier to establish which gloves are suitable for such an application.

Bone cement application:

- Clinical studies show the need to maintain strictly aseptic surgical procedures. Any deep infection of a surgical wound is a serious risk and will affect the successful outcome of the technique. Deep wound infection is a serious post-operative complication and may require total removal of the embedded cement. Deep wound infec-

tion may be latent and not manifest itself even for several years post-operatively.

- Patients should be carefully monitored for any change in blood pressure during and immediately following the application of bone cement. Adverse patient reactions affecting the cardiovascular system have been associated with the use of bone cements, these include: hypotension, hypoxemia, cardiac arrhythmia, bronchospasm, cardiac arrest, myocardial infarction, pulmonary embolism, cerebrovascular accident and possible death. Hypotensive reactions have occurred between 10 and 165 seconds following application of bone cement; they have lasted from 30 seconds to 5 or more minutes. Some have progressed to cardiac arrest. In addition, the over-pressurization of the bone cement should be avoided during the insertion of the bone cement and the implant in order to minimize the occurrence of pulmonary embolism.
- The preparation of the bone marrow cavity results in marrow contents entering the blood stream. Prior to the application of bone cement to the bone, the cavity should be thoroughly cleaned by brushing and washing (lavage) to remove fat, marrow and other debris. The cavity should be kept as dry as possible to prevent blood and debris becoming mixed with the cement. Thorough cleaning of the bone reduces the risk of marrow content being forced into the vascular system during the insertion of bone cement and subsequent pressurization. The expulsion of bone marrow has been associated with the occurrence of pulmonary embolisms, and this risk has been found to be increased in patients with highly osteoporotic bone and patients diagnosed with femoral neck fracture. Reaming of the marrow cavity can have similar effects on mean arterial pressure as the introduction of the bone cement. Marrow cavities should be vented when the cement is introduced digitally.
- Consideration should be given to the use of acrylic bone cement in patients diagnosed with femoral neck fracture, as some published literature has indicated there is a potential for increased mortality compared with uncemented techniques.
- Inadequate fixation or unanticipated post-operative events may affect the cement-bone interface and lead to micromotion of the cement against bone surfaces with which the cement is in contact. A fibrous tissue layer may develop between the cement and the bone. Long term follow-up is advised for all patients on a regular scheduled basis.
- Prosthesis to be implanted must be compatible with the use of bone cement.
- The completion of cement polymerization occurs in the patient and is an exothermic reaction with considerable liberation of heat. The long term effects of the heat produced in situ have not yet been established.

Incompatibilities

- Aqueous solutions (e.g. ones containing antibiotics) must not be added to the bone cement because they have a considerable detrimental effect on the physical and mechanical properties of the cement.

Precautions during pregnancy, breast-feeding and in children.

- The safety and effectiveness of **G1 40** bone cements in pregnant women or in children has not yet been established.

Adverse events

Serious adverse events, some with fatal outcome, associated with the use of acrylic bone cements include:

- Myocardial infarction.
- Cardiac arrest.
- Cerebrovascular accident.
- Pulmonary embolism.
- Anaphylaxis.

The most frequent adverse reactions reported with acrylic bone cements are:

- Transitory fall in blood pressure.
- Elevated serum gamma-glutamyl-transpeptidase (GGTP) up to 10 days post-operation.
- Thrombophlebitis.
- Hemorrhage and hematoma.
- Pain and/or loss of function.
- Loosening or displacement of the prosthesis.
- Superficial or deep wound infection.
- Trochanteric bursitis.
- Short-term cardiac conduction irregularities.
- Heterotopic new bone formation.
- Trochanteric separation.

Other potential adverse events reported for bone cements include:

- Hypoxemia.
- Cardiac arrhythmia.
- Bronchospasm.
- Adverse tissue reaction.
- Pyrexia due to allergy to the bone cement.
- Hematuria.
- Dysuria.
- Bladder fistula.
- Local neuropathy.
- Local vascular erosion and occlusion.
- Transitory worsening of pain due to heat released during polymerization.
- Delayed sciatic nerve entrapment due to extrusion of the bone cement beyond the region of its intended application.
- Intestinal obstruction because of adhesions and stricture of the ileum due to the heat released during cement polymerization.

Dose

Choose the dose considering the intervention to be performed. As a preventive measure, it is advisable to keep at hand at least a second package of **G1 40** bone cement to be able to prepare an extra quantity of bone cement, if required during the procedure.

The maximum recommended dose of **G1 40** bone cement is of 3 packs.

Directions for use

Cement Preparation

- The protective outer foil pouch, the outer peelable pouch of the powder component and the blister pack enclosing the ampoule of liquid component, should be opened by a circulating nurse. The inner bag (or pouch) containing the powder component and the sterile ampoule containing the liquid component are aseptically transferred into the sterile operative area.
- The sterile bag (or pouch) containing the powder component is opened with sterile scissors and the entire contents are emptied into a suitable clean, dry, sterile mixing vessel made from an inert material (such as glass, ceramic, stainless steel, or non-reactive plastics). The sterile ampoule containing the liquid component is opened and the entire contents are emptied evenly onto the powder in the mixing vessel.
- A standard dose of bone cement is prepared by mixing the entire liquid contents of the ampoule with the entire contents of the powder bag. The amount of mixed cement required for clinical use is determined by the surgeon in each individual case.

Mixing and Digital Application

- **G1 40** can be applied digitally. The surgeon must use their clinical judgement to decide when the cement is of a suitable viscosity to allow the surgical procedure to continue.
- Prior to cement application, it is recommended that a cement restrictor is always used during cementation of the femur and that this is introduced at the required depth.
- The cement is mixed thoroughly but carefully to minimize the entrapment of air. Once a dough is formed the surgeon should wait until the cement no longer adheres to the glove. The cement can then be taken into gloved hands and kneaded thoroughly. It is vital that premature insertion of cement is avoided as this may lead to a drop in the patient's blood pressure. To avoid this, the appearance of the cement should be observed to ensure the surface has become dull as opposed to shiny. Also cement should not adhere excessively to the surgeon's gloves. The time of cement application and prosthesis insertion is at the discretion of the surgeon and will depend upon the surgical procedure used.
- Implant insertion should be carried out at a time appropriate for the bone/joint and prosthesis design concerned. In general, implant insertion should be delayed until the cement has developed a sufficient degree of viscosity to resist excessive displacement by the implant. However, implant insertion should not be delayed such that there is a risk that the procedure cannot be completed due to cement hardening.
- Following introduction the implant must be firmly held in position to avoid movement and pressurization must be maintained until the cement finally hardens. Excess bone cement must be removed before the cement has completely hardened.

The handling characteristics and setting times are affected by ambient temperature. Please refer to the end of the instruction leaflet for guidance charts (Note: the usage charts were generated under controlled laboratory conditions).

Storage

- Store the sealed outer pack at a maximum temperature of 77°F (25°C) and protect it from light to prevent premature polymerization of the liquid monomer component.
- Store at the intended operating room temperature for a minimum 24 hrs before use.
- Do not use the product after the expiration date.

Disposal

Dispose the remaining of **G1 40** bone cement and the content of partially utilized, expired, not usable or damaged packages in an authorized waste facility; any leftover liquid component should be evaporated under a well-ventilated hood or absorbed by an inert material and transferred to a suitable container for disposal.

G3 40 Radiopaque bone cement

Device Description

G3 40 bone cement is a self-curing, radiopaque, polymethylacrylate (PMMA) based cement which is used for securing a metal or polymeric prosthesis to living bone in arthroplasty procedures. The hardened bone cement secures the fixation of the grafted artificial joint improving the transfer of forces at the interface of the implant-bone.

The bone cement is supplied as a two-component system, consisting of separate, sterile liquid and powder components, which are mixed together at the point of use to produce the cement.

The liquid component is sterilized by membrane filtration and aseptically filled into a sterile glass ampoule. The ampoule is contained within a sealed blister pack, which is sterilized using ethylene oxide. The powder component is contained in a Tyveck/PE pouch, within a peelable pouch and is sterilized by ethylene oxide. The sterile powder component is supplied within an outer, protective, non-sterile, laminated foil pouch.

G3 40 bone cement is a low viscosity cement, primarily intended for syringe application.

Pack size	
Powder weight (g)	34,40
Liquid volume (ml)	20,00

Composition

Composition of the powder component (% w/w)	
Polymethyl Methacrylate	88,8
Benzoyl Peroxide	10,0
Barium Sulfate	1,2

Composition of the liquid component (% w/w)	
Methyl Metacrylate	97,50
N,N dimethyl-p-toluidine	2,5
Methy Ether of Hydroquinone	50 ppm

Intended use

G3 40 radiopaque bone cements is intended for use in arthroplasty procedures of the hip, knee, ankle, shoulder and other joints for the fixation of polymer or metallic prosthetic implants to living bone.

Contraindications

The use of **G3 40** is contraindicated:

- when a patient is known to be hypersensitive to the constituents of bone cement or to the contrast medium (barium sulfate);
- when the local or systemic infections not completely resolved;
- where there is no possibility of bone regeneration;
- during pregnancy or breast-feeding: The safety and effectiveness of the radiopaque bone cements in pregnant women or in children has not yet been established.

Precautions/Warnings

Bone cement preparation:

- **G3 40** bone cement is supplied sterile for single use only. Do not re-use. Re-sterilization of any components of the cement must not be attempted.
- Ensure the inner packages and components are undamaged. Powder should be consistent (no agglomerations) and not yellow or brown in color. The contents within the vial should appear as a low viscosity liquid. Do not use the liquid monomer if it shows any sign of thickening or premature polymerization. Always check the condition of the liquid monomer before performing the procedure. If the powder has a yellowish or brownish color or if the liquid is syrupy, do not use product. This indicates the product has not been stored properly.
- A dose is prepared by mixing the entire contents of the bag of cement powder with all the monomer liquid of an ampoule. The quantity of cement dough required depends on the specific surgical intervention and on the technique being used. At least one additional dose of **G3 40** should be available before commencing the operation.
- The protective outer foil pouch and the outer peelable pouch are non-sterile and must not be transferred into the sterile field
- The surgeon should, by specific training and experience, be thoroughly familiar with the properties, handling characteristics and application of bone cements. The user is advised to practice the entire procedure of mixing, handling and introducing **G3 40** before using it for the first time. Detailed knowledge is necessary even if mixing systems and syringes are being used for application of the cement. Because the handling and curing characteristics of these cements vary with temperature and mixing technique, they are best determined by the surgeon's actual experience.
- Handling characteristics and setting time of bone cements are affected by:
 - o temperature
 - o humidity
 - o vacuum mixing
- Refer to clinical timing and usage charts at the end of this leaflet for effect of temperature on handling and setting times.

- Manual handling of the cement and body temperature will reduce the final setting time.
- Variations in humidity will affect the cement handling characteristics and setting time.
- Variations in setting time over the cement's shelf life can be minimized by storing the cement under the recommended conditions
- The surgeon should be familiar with the effects of a particular mixing system on handling characteristics and setting time of the cement before use. Vacuum mixing can noticeably accelerate the setting time of the product.

Handling of the liquid monomer:

- As the liquid monomer is highly volatile and flammable, the operating room should be adequately ventilated to eliminate as much monomer vapor as possible. Ignition of monomer fumes caused by use of electrocautery devices in surgical sites near freshly implanted bone cements has been reported.
- Caution should be exercised during the mixing of the two components to prevent excessive exposure to the concentrated vapors of the monomer, which may produce: irritation of the respiratory tract, eyes, and possibly the liver; contact dermatitis. If the liquid component comes into contact with the eyes, wash with copious amounts of water.
- Concentrated vapors of the liquid component may have an adverse reaction with contact lenses. Since soft contact lenses are permeable to liquids and gases, they should be worn in the operating theater if methyl methacrylate is being used.
- Methyl methacrylate has been demonstrated to cause hypersensitivity in susceptible persons, which may result in an anaphylactic response.
- The liquid component of bone cement is a powerful lipid solvent. This liquid component should not be allowed to come into contact with rubber or latex surgical gloves. Wearing of a second pair of gloves and strict adherence to the mixing instructions may diminish the possibility of hypersensitivity reactions. Gloves made of PVP (polyethylene, ethylene vinylalcohol copolymer, polyethylene) and Viton ®/butyl gloves have proved to provide good protection over a lengthy period. For safety's sake it is recommended that two pairs of gloves be worn over one another, e.g. one polyethylene surgical glove over an inner pair of standard latex surgical gloves. The use of latex or polystyrene-butadiene gloves on their own is inadequate. Please make enquiries with your supplier to establish which gloves are suitable for such an application.

Bone cement application:

- Clinical studies show the need to maintain strictly aseptic surgical procedures. Any deep infection of a surgical wound is a serious risk and will affect the successful outcome of the technique. Deep wound infection is a serious post-operative complication and may require total removal of the embedded cement. Deep wound infection

may be latent and not manifest itself even for several years post-operatively.

- Patients should be carefully monitored for any change in blood pressure during and immediately following the application of bone cement. Adverse patient reactions affecting the cardiovascular system have been associated with the use of bone cements, these include: hypotension, hypoxemia, cardiac arrhythmia, bronchospasm, cardiac arrest, myocardial infarction, pulmonary embolism, cerebrovascular accident and possible death. Hypotensive reactions have occurred between 10 and 165 seconds following application of bone cement; they have lasted from 30 seconds to 5 or more minutes. Some have progressed to cardiac arrest. In addition, the over-pressurization of the bone cement should be avoided during the insertion of the bone cement and the implant in order to minimize the occurrence of pulmonary embolism.
- The preparation of the bone marrow cavity results in marrow contents entering the blood stream. Prior to the application of bone cement to the bone, the cavity should be thoroughly cleaned by brushing and washing (lavage) to remove fat, marrow and other debris. The cavity should be kept as dry as possible to prevent blood and debris becoming mixed with the cement. Thorough cleaning of the bone reduces the risk of marrow content being forced into the vascular system during the insertion of bone cement and subsequent pressurization. The expulsion of bone marrow has been associated with the occurrence of pulmonary embolisms, and this risk has been found to be increased in patients with highly osteoporotic bone and patients diagnosed with femoral neck fracture. Reaming of the marrow cavity can have similar effects on mean arterial pressure as the introduction of the bone cement. Marrow cavities should be vented when the cement is introduced digitally.
- Consideration should be given to the use of acrylic bone cement in patients diagnosed with femoral neck fracture, as some published literature has indicated there is a potential for increased mortality compared with uncemented techniques.
- Inadequate fixation or unanticipated post-operative events may affect the cement-bone interface and lead to micromotion of the cement against bone surfaces with which the cement is in contact. A fibrous tissue layer may develop between the cement and the bone. Long term follow-up is advised for all patients on a regular scheduled basis.
- Prosthesis to be implanted must be compatible with the use of bone cement.
- The completion of cement polymerization occurs in the patient and is an exothermic reaction with considerable liberation of heat. The long term effects of the heat produced in situ have not yet been established.

Incompatibilities

- Aqueous solutions (e.g. ones containing antibiotics) must not be added to the bone cement because they have a considerable detrimental effect on the physical and mechanical properties of the cement.

Precautions during pregnancy, breast-feeding and in children.

- The safety and effectiveness of **G3 40** bone cements in pregnant women or in children has not yet been established.

Adverse events

Serious adverse events, some with fatal outcome, associated with the use of acrylic bone cements include:

- Myocardial infarction.
- Cardiac arrest.
- Cerebrovascular accident.
- Pulmonary embolism.
- Anaphylaxis.

The most frequent adverse reactions reported with acrylic bone cements are:

- Transitory fall in blood pressure.
- Elevated serum gamma-glutamyl-transpeptidase (GGTP) up to 10 days post-operation.
- Thrombophlebitis.
- Hemorrhage and hematoma.
- Pain and/or loss of function.
- Loosening or displacement of the prosthesis.
- Superficial or deep wound infection.
- Trochanteric bursitis.
- Short-term cardiac conduction irregularities.
- Heterotopic new bone formation.
- Trochanteric separation.

Other potential adverse events reported for bone cements include:

- Hypoxemia.
- Cardiac arrhythmia.
- Bronchospasm.
- Adverse tissue reaction.
- Pyrexia due to allergy to the bone cement.
- Hematuria.
- Dysuria.
- Bladder fistula.
- Local neuropathy.
- Local vascular erosion and occlusion.
- Transitory worsening of pain due to heat released during polymerization.
- Delayed sciatic nerve entrapment due to extrusion of the bone cement beyond the region of its intended application.
- Intestinal obstruction because of adhesions and stricture of the ileum due to the heat released during cement polymerization.

Dose

Choose the dose considering the intervention to be performed. As a preventive measure, it is advisable to keep at hand at least a second package of **G3 40** bone cement to be able to prepare an extra quantity of bone cement, if required during the procedure.

The maximum recommended dose of **G3 40** bone cement is n° 3 packs

Directions for use

Cement Preparation

- The protective outer foil pouch, the outer peelable pouch of the powder component and the blister pack enclosing the ampoule of liquid component, should be opened by a circulating nurse. The inner bag (or pouch) containing the powder component and the sterile ampoule containing the liquid component are aseptically transferred into the sterile operative area.
- The sterile bag (or pouch) containing the powder component is opened with sterile scissors and the entire contents are emptied into a suitable clean, dry, sterile mixing vessel made from an inert material (such as glass, ceramic, stainless steel, or non-reactive plastics). The sterile ampoule containing the liquid component is opened and the entire contents are emptied evenly onto the powder in the mixing vessel.
- A standard dose of bone cement is prepared by mixing the entire liquid contents of the ampoule with the entire contents of the powder bag. The amount of mixed cement required for clinical use is determined by the surgeon in each individual case.

Mixing and Syringe application

- **G3 40** bone cements may be applied using a suitable cement gun and syringe.
 - The bone cement is prepared and mixed as described previously by adding all of the liquid component to all the powder component. The cement is then transferred into a suitable cement gun cartridge. The surgeon should use their experience to judge when the cement has reached an appropriate viscosity to be extruded. This will not occur until after the cement has formed a dough. A small amount of cement should be extruded from the syringe and visually assessed to ensure that the surface of the cement appears dull and excessive flow under gravity has ceased.
 - Prior to extrusion, it is recommended that a cement restrictor be inserted, at the required depth into, the prepared bone cavity. Introduction of bone cement into the prepared cavity should be carried out in a retrograde fashion. Once the cavity is filled it is strongly advised that adequate pressurization is applied and maintained up to the point of hardening. Implant insertion should be carried out at a time appropriate for the bone/joint and prosthesis design concerned. In general, implant insertion should be delayed until the cement has developed a sufficient degree of viscosity to resist excessive displacement by the implant. However, implant insertion should not be delayed such that there is a risk that the procedure cannot be completed due to
 - cement hardening.
- Following introduction the implant must be firmly held in position to avoid movement and pressurization must be maintained until the cement finally hardens. Excess bone cement must be removed before the cement has completely hardened

The handling characteristics and setting times are affected by ambient temperature. Please refer to the end of the instruction leaflet for guidance charts (Note: the usage charts were generated under controlled laboratory conditions).

Storage

- Store the sealed outer pack at a maximum temperature of 77°F (25°C) and protect it from light to prevent premature polymerization of the liquid monomer component.
- Store at the intended operating room temperature for a minimum 24 hrs before use.
- Do not use the product after the expiration date.

Disposal

Dispose the remaining of **G3 40** bone cement and the content of partially utilized, expired, not usable or damaged packages in an authorized waste facility; any leftover liquid component should be evaporated under a well-ventilated hood or absorbed by an inert material and transferred to a suitable container for disposal.

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