



IMPLANTING POSSIBILITIES. IMPROVING LIVES.


About Company

ShalbyMedTech medical devices ltd is a leading Orthopedic implant company with a diverse product portfolio and rapidly expanding global footprints. As the leading manufacturer and supplier of medical devices, ShalbyMedTech medical devices ltd has attained recognition globally by continuing its focus on creating precise and safe orthopedic products. Our range of orthopedic implants and prostheses helps surgeons offer effective orthopedic treatments, shorter hospitalizations, and better overall quality of life for patients.

We believe in the continuous development of our product portfolio and regularly launch new and innovative orthopedic products to meet the changing landscape of orthopedic technologies and the needs of surgeons.


Femoral Metal Head

Unique Locking Mechanism



Metal Head			
Cat. No.	Specification	Product Description	Diameter
MNH1022S	22S	Product Name: Metal Head JX 4301 Material: CoCrMo	φ 22
MNH1022M	22M		φ 22
MNH1022L	22L		φ 22
MNH1028S	28S		φ 28
MNH1028M	28M		φ 28
MNH1028L	28L		φ 28
MNH1028XL	28XL		φ 28
MNH1028XXL	28XXL		φ 28
MNH1032S	32S		φ 32
MNH1032M	32M		φ 32
MNH1032L	32L		φ 32
MNH1032XL	32XL		φ 32
MNH1032XXL	32XXL		φ 32
MNH1036S	36S		φ 36
MNH1036M	36M		φ 36
MNH1036L	36L		φ 36
MNH1036XL	36XL		φ 36

CEMSET™ CEMENTED FEMORAL STEM



CEMSET™ CEMENTED FEMORAL STEM

Cat. No.	Specification	Product Description	Neck Shaft Angle	Offset (B)	Neck Length (C)	Stem length (L)	Distal Dia.
MSC40505	1#	Product Name: RCH Cemented Femoral Stem JX 1401H Material: CoCrMo Surface: High Polished	130°	36.4	35	120	7
MSC40506	2#			38.4	37	125	7
MSC40507	3#			38.4	37	130	7.5
MSC40508	4#			40.5	39	140	8
MSC40509	5#			40.5	39	145	8.5
MSC40510	6#			42.4	41	150	9
MSC40511	7#			42.4	41	155	10

CEMSET™ LONG CEMENTED FEMORAL STEM

Cat. No.	Specification	Product Description	Neck Shaft Angle	Offset (B)	Neck Length (C)	Stem Length (L)	Distal Dia. (mm)
MSC40507L	3#	Product Name: RCH Long Cemented Femoral Stem JX 1403H Material: CoCrMo Surface: High Polished	130°	38.4	37	180	8
MSC40508L	4#			40.5	39	190	8
MSC40509L	5#			40.5	39	200	8.5
MSC40510L	6#			42.4	41	210	9
MSC40511L	7#			42.4	41	220	9

CEMCTM CEMENTED ACETABULAR CUP

Cat. No.	Specification	Product Description	Inner Dia. (mm)
MC6040	40/28	Product Name: Cemented Acetabular JX 21201 Material: UHMWPE (Ultrahigh molecular weight polyethylene)	28
MC6042	42/28		28
MC6044	44/28		28
MC6046	46/28		28
MC6048	48/28		28
MC6050	50/28		28
MC6052	52/28		28
MC6054	54/28		28
MC6056	56/28		28
MC6058	58/28		28



SHALBY
MedTech

Shalby Hospitals, Opp: Karnawati Club, S. G. Highway, Near
Prahladnagar Garden, NA, AHMEDABAD, Ahmedabad-380015, Gujarat, India
Email : infomedtech@shalby.org | Phone: +91 6357007339