



SHALBY
MedTech

IMPLANTING POSSIBILITIES. IMPROVING LIVES.

REVSET™

REVSET™

CEMENTLESS FEMORAL STEM

- Forged titanium alloy, and monobloc conical stem.
- Designed for secure fixation by means of longitudinal ribs and 2 taper and allow for intraoperative anteversion.
- Eight longitudinal ribs provide rotational stability and anti-rotation effect. Varying implant lengths are available (190–265 mm).
- Simple and precise tools.
- Clinically proven design with over 20 years of clinical experience.
- The revision stent is recommended as the first choice by the Swedish Arthroplasty Register.



REVSET™ CEMENTLESS FEMORAL STEM



Cat. No.	Specification	Product Description	Neck Shaft Angle	Offset (B)	Length (C)	Length (L)	Dia. (D)	Distal Dia. (F)
M42314-190L	1#/190	Product Name: RWH Cementless Femoral Stem JX F1103A Material: Ti6Al4V Surface: Ti Coating	135°	42	59	190	14	10.4
M42315-190L	2#/190			42	59		15	11.4
M42316-190L	3#/190			42	59		16	12.4
M42317-190L	4#/190			42	59		17	13.4
M42318-190L	5#/190			44	62		18	14.4
M42319-190L	6#/190			44	62		19	15.4
M42320-190L	7#/190			44	62		20	16.4
M42321-190L	8#/190			42	59		12	8.4
M42322-190L	9#/190			42	59		13	9.4
M42314-225L	1#/225			42	59	225	14	10.4
M42315-225L	2#/225			42	59		15	11.4
M42316-225L	3#/225			42	59		16	12.4
M42317-225L	4#/225			42	59		17	13.4
M42318-225L	5#/225			44	62		18	14.4
M42319-225L	6#/225			44	62		19	15.4
M42320-225L	7#/225			44	62		20	16.4
M42321-225L	8#/225			44	62		21	17.4
M42322-225L	9#/225			46	64		22	18.4
M42314-265L	1#/265			42	59	265	14	10.4
M42315-265L	2#/265			42	59		15	11.4
M42316-265L	3#/265			42	59		16	12.4
M42317-265L	4#/265			42	59		17	13.4
M42318-265L	5#/265			44	62		18	14.4
M42319-265L	6#/265			44	62		19	15.4
M42320-265L	7#/265			44	62		20	16.4
M42321-265L	8#/265			44	62		21	17.4
M42322-265L	9#/265			46	65		22	18.4
M42323-265L	10#/265			46	65		23	19.4

ABOUT COMPANY

ShalbyMedTech medical devices Ltd is a leading Orthopedic implant company with a diverse product portfolio and rapidly expanding global footprints. As the leading manufacturer and supplier of medical devices, ShalbyMedTech medical devices Ltd has attained recognition globally by continuing its focus on creating precise and safe orthopedic products. Our range of orthopedic implants and prostheses helps surgeons offer effective orthopedic treatments, shorter hospitalizations, and better overall quality of life for patients.

We believe in the continuous development of our product portfolio and regularly launch new and innovative orthopedic products to meet the changing landscape of orthopedic technologies and the needs of surgeons.



Address: Shalby Hospitals, Opp: Karnawati Club, S. G. Highway, Near
Prahladnagar Garden, NA, AHMEDABAD, Ahmedabad-380015, Gujarat, India

Email: infomedtech@shalby.org | **Phone:** +91 6357007339